

Enhancing Sequence Coverage of Integral Membrane Proteins Using Proteomic Approaches: Application to the G-protein Coupled Receptor for Neurotensin



Jenny TC Ho & Sonja Hess

Mass Spectrometry and Proteomics Facility, NIDDK, NIH, DHHS, Bethesda, MD

Overview

- Apply mass spectrometry based proteomic approaches to help understand the mechanisms involved in the recognition, binding and activation of the G-protein coupled receptor (GPCR) for neurotensin (NTR1).
- To obtain complete sequence coverage of rat NTR1 to monitor conformational changes, characterize mutations in protein sequence and post-translational modifications.

Introduction: GPCRs

- GPCRs are involved in signal recognition and cell communication. Over 1,000 different GPCRs have been identified and many are implicated as therapeutic drug targets.
- Neurotensin is involved in modulatory functions for example in the heart, digestive system and in the central nervous system.
- When the target ligand (for example neurotensin) binds to its GPCR, the receptor undergoes a conformational change, this in turn activates G-proteins and consequently activates or inhibits effector molecules.




Figure 1: Ribbon drawing of GPCR rhodopsin²

- GPCRs are the largest family of receptors. They all have a similar structural framework; seven helices that traverse the lipid membrane linked by intra- and extra-cellular loops¹ (figure 1).
- Analysis of these receptors by MS based proteomic methods are more challenging compared to soluble proteins due to their hydrophobicity.
- These proteins are prone to aggregation and require detergents for solubilization, which can have deleterious effects on LC separation and mass spectral analysis.

References

1. T.P. Iismaa *et al.*, G Protein-Coupled Receptors, R.G. Landes Company, 1995
2. http://blanco.biomol.uci.edu/Membrane_Proteins_xtal.html
3. J.F. White *et al.*, FEBS Letters, 2004, 564, 289-293.

Acknowledgments

Many thanks to R. Grishammer (NIDDK, NIH) for providing NTR1 and to G. Abdoulaeva, CBER, FDA for synthesizing transmembrane peptide.

Experimental




Figure 2: Flow diagram outlining the procedures used to obtain peptide mass map of rat NTR1.

- Rat NTR1 was expressed and purified³. Figure 2 shows the procedures used to obtain peptide mass map of NTR1. LC-MS/MS was performed using a NanoAquity LC system (C18 Asymmetry column, ACN containing 0.2% formic acid gradient) and Q-ToF II MS instrument (Water/Micromass).
- Transmembrane peptide mimic (TMpep 1-66) was synthesized. The crude material was solubilized in 50% formic acid and purified on a Zorbax C3 column (Agilent) employing an ACN gradient containing 0.2% formic acid. ES analysis of transmembrane peptides was performed by direct infusion using a Nanomate (Advion) into an LTQ-FTICR (Thermo).

Results I – Optimization of methods using transmembrane peptide mimics

- A bottom-up proteomics approach was employed to obtain complete sequence coverage of NTR1.
- Insufficient sequence coverage of NTR1 was obtained when samples prepared following procedures generally used in proteomic studies. For example, only 15% and 8% of rat NTR1 was observed when digested with trypsin in 8M urea and with proteinase K respectively (figure 3).




Figure 3: Amino acid sequence of rat NTR1. Yellow and red circles represent residues observed by LC-MS/MS after trypsin digestion in urea and proteinase K digestion, respectively.

- Secondary structure, solubility, peptide length and hydrophobicity may explain why peptides that make up the transmembrane (TM) regions are rarely observed in proteomic analyses. To evaluate these parameters, a peptide that resembles TM regions was synthesized. The AA sequence of the TM peptide mimic (TMpep1-66) and ES spectrum is shown in figure 4. This sequence corresponds to TM I and II of rat NTR1.




Figure 4: Amino acid sequence of transmembrane peptide mimic (TMpep 1-66). ES spectrum of purified TMpep 1-66, spray solvent 80/20 ACN/H₂O, 0.2% formic acid

- Ammonium bicarbonate with the addition of organic solvents eg methanol, acetonitrile or acetone, were investigated to solubilize TMpep 1-66. Organic solvents were chosen that were also compatible with enzymatic digestion, RPLC and MS analysis.
- Figure 5 shows selected ion chromatograms of peptides generated from a trypsin digest of TMpep 1-66 in the presence of 60% methanol. Successful RPLC separation of hydrophobic peptides that make up TM domains was achieved. Peak tailing was more prominent with increasing peptide hydrophobicity.




Figure 5: LC-MS analysis of TMpep 1-66 tryptic digest. Selected ion chromatogram of TMpep 1-66³⁺, TMpep 1-31³⁺ (1 missed cleavage), TMpep 5-31³⁺ and TMpep 34-66³⁺. Dotted line indicates % of ACN gradient containing 0.2% formic acid.




Figure 6: ES MS/MS spectrum of tryptic peptides generated from the digestion of TMpep 1-66 (see figure 4). (a) MS/MS spectrum of TMpep 5-31³⁺ VLVTAIYLALFVVGTVGNSVTAFTLAR (makes up TM I of rat NTR1). (b) MS/MS spectrum TMpep 34-66³⁺ SLQSLQSTVHYHLGSLALSDLILLAMPVELY (makes up TM II of rat NTR1).

- Tryptic peptides; TMpep 5-31 and TMpep 34-66 were subjected to CAD in a Q-ToF and LTQ-FT instrument. Figure 6 shows MS/MS spectra of M+3H³⁺ TM peptides acquired using a LTQ-FTICR.

Results II – Mapping NTR1




Figure 7: Amino acid sequence of rat NTR1. Blue, yellow and green circles represent residues observed by LC-MS/MS after trypsin+chymotrypsin digestion in Rapigest, CNBr followed by trypsin digestion and in both digestion methods, respectively.

- Methods optimized using a synthetic TM peptide mimic were applied to recombinant rat NTR1. Figure 7 shows the AA sequence of rat NTR1, where the colored circles indicate peptides observed by LC-MS/MS either after digestion with CNBr followed by trypsin or trypsin+chymotrypsin digestion in 5% acid labile surfactant (Rapigest, Waters).
- Peptides corresponding to 5 out of the 7 TM regions have been observed, as well as the intra- and extra-cellular loops.

Conclusions

- The choice of solvent used to prepare integral membrane protein for mass spectral analysis is important. It should (i) solubilize the protein (ii) be compatible with enzymatic/chemical digestion and (iii) be compatible with down stream analysis; LC separation and mass spectral analysis.
- To improve our understanding as to why transmembrane peptides are often under-represented in proteomic analyses, a synthetic peptide was used to optimize sample preparation (i.e test various solubilizing reagents and digestion buffers) and also LC-MS/MS analyses.
- Precipitation of NTR1 reduced concentration of detergents. High % of acid was more effective in solubilizing the precipitated protein compared to organic solvents.
- Using these conditions we are now able to increase sequence coverage from 24% to 60% of the rat NTR1 protein when using multiple enzymes (figure 7).